

## **QuiKResponse - Low-Level Mercury, waste water (WasteWatR TM)**

Art. ID                      ERA-931QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 5 mL flame-sealed ampule yields up to 4 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Mercury (Hg)	[7439-97-6]	1 - 100	ng/L		