

QuiKResponse - Low-Level Mercury, waste water (WasteWatR TM)

Art. ID ERA-931QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 5 mL flame-sealed ampule yields up to 4 liters after dilution.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Mercury (Hg) | [7439-97-6] | 1 - 100 | ng/L | | |