

## **QuiKResponse - pH, waste water (WasteWatR TM)**

Art. ID ERA-977QR  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 250 mL whole-volume bottle is ready to analyze. Use with electrometric methods.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	pH		5,0 - 10,0	S.U.		