

## **QuiKResponse - pH, waste water (WasteWatR TM)**

Art. ID ERA-977QR  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 250 mL whole-volume bottle is ready to analyze. Use with electrometric methods.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|------------|---------------------|------|--------|--------|
|                  | pH                |            | 5,0 - 10,0          | S.U. |        |        |