

## **QuiKResponse - Ignitability/Flash Point, priority pollutant (PriorityPollutnT TM)**

Art. ID                      ERA-979QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One standard packaged in three 30 mL bottles. Use to measure ignitability.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Flash point		100 - 200	°F-closed cup		