

QuiKResponse - Ignitability/Flash Point, priority pollutant (PriorityPollutnT TM)

Art. ID ERA-979QR
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One standard packaged in three 30 mL bottles. Use to measure ignitability.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Flash point		100 - 200	°F-closed cup		