

QuiKResponse - Corrosivity, potable water (PotableWatR TM)

Art. ID	ERA-980QR
Unit	each
Deliverydetails	No Dangerous Good /not restricted

Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 500 mL whole-volume bottle is ready to analyze for corrosivity, calcium carbonate saturation and Langelier saturation index.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Corrosivity		- 4 to +4	SI units		