

## **QuiKResponse - Cyanide, potable water (PotableWatR TM)**

Art. ID                      ERA-983QR  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Free cyanide		0,1 - 0,5	mg/L		