

## **QuiKResponse - Hexavalent Chromium, waste water (WasteWatR TM)**

Art. ID ERA-984QR  
Unit each  
Deliverydetails No Dangerous Good /not restricted

### Description

On-demand QuiKResponse PTs - the quickest turn in the industry and available 52 weeks a year! One 15 mL screw-cap vial yields up to 2 liters after dilution. Use with IC or colorimetric methods.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Hexavalent chromium	[18540-29-9]	45 - 880	µg/L		