

Zinc - DSC calibration standard

| | |
|-----------------|-----------------------------------|
| Art. ID | LGC2611 |
| Unit | 500 mg |
| Deliverydetails | No Dangerous Good /not restricted |

Description

This certified reference material is intended for use in the calibration of differential scanning calorimeters and similar instruments. Certified values:

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---------------------|------------|---------------------|--------|--------|--------|
| | Enthalpy of Fusion | | 7,103 | kJ/mol | | |
| | Melting temperature | | 419,53 | °C | | |