

## **Whisky - Congeners**

Art. ID LGC5100  
Unit 10 ml  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1170 III E1

### **Description**

Whisky sourced from a commercial source is available in 10mL unit in amber glass vials using sealed septa with crimped caps. Indicative values for acetaldehyde, butan-1-ol, furfural, ethyl acetate. Certified values:

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Methanol	[67-56-1]	5,2	g/100 L		
	1-Propanol	[71-23-8]	57	g/100 L		
	2-Methyl-1-propanol	[78-83-1]	58,8	g/100 L		
	2-Methyl-1-butanol	[137-32-6]	21,38	g/100 L		
	3-Methyl-1-butanol	[123-51-3]	58,2	g/100 L		
	Butanol	[71-36-3]	0,48	g/100 L		