

Waste water - Trace metals

Art. ID SPS-WW2
Unit 6 x 50 mL
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	10000	ng/mL		
	Arsenic (As)	[7440-38-2]	500	ng/mL		
	Cadmium (Cd)	[7440-43-9]	100	ng/mL		
	Cobalt (Co)	[7440-48-4]	300	ng/mL		
	Chromium (Cr)	[7440-47-3]	1000	ng/mL		
	Copper (Cu)	[7440-50-8]	2000	ng/mL		
	Iron (Fe)	[7439-89-6]	5000	ng/mL		
	Manganese (Mn)	[7439-96-5]	2000	ng/mL		
	Nickel (Ni)	[7440-02-0]	5000	ng/mL		
	Phosphorus (P)	[7723-14-0]	5000	ng/mL		
	Lead (Pb)	[7439-92-1]	500	ng/mL		
	Vanadium (V)	[7440-62-2]	500	ng/mL		
	Zinc (Zn)	[7440-66-6]	3000	ng/mL		