

Aluminium base, Si, RM, Chips 100 g

Art. ID BAS-BCS-CRM 380/1
Unit 100 g
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Copper (Cu) | [7440-50-8] | 0,91 | % | | |
| | Iron (Fe) | [7439-89-6] | 1,24 | % | | |
| | Magnesium (Mg) | [7439-95-4] | 0,24 | % | | |
| | Manganese (Mn) | [7439-96-5] | 0,094 | % | | |
| | Nickel (Ni) | [7440-02-0] | 0,94 | % | | |
| | Silicon (Si) | [7440-21-3] | 1,93 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,024 | % | | |
| | Zinc (Zn) | [7440-66-6] | 0,025 | % | | |