

### **Copper base, Gun Metal, CRM, Chips 100 g**

Art. ID                      BAS-BCS-CRM 207/2  
Unit                         100 g (chips)  
Deliverydetails            No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,013	%		
	Arsenic (As)	[7440-38-2]	0,066	%		
	Bismuth (Bi)	[7440-69-9]	0,04	%		
	Copper (Cu)	[7440-50-8]	87,35	%		
	Iron (Fe)	[7439-89-6]	0,029	%		
	Nickel (Ni)	[7440-02-0]	0,28	%		
	Lead (Pb)	[7439-92-1]	0,7	%		
	Antimony (Sb)	[7440-36-0]	0,1	%		
	Silicon (Si)	[7440-21-3]	0,016	%		
	Tin (Sn)	[7440-31-5]	9,74	%		
	Zinc (Zn)	[7440-66-6]	1,6	%		