

Iron base, Cast Iron, RM, Block 40 mm x 37 mm x 10 mm

Art. ID BAS-LARM 2
Unit block
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,066	%		
	Arsenic (As)	[7440-38-2]	0,044	%		
	Carbon (C)	[7440-44-0]	~3	%		
	Cerium (Ce)	[7440-45-1]	0,008	%		
	Chromium (Cr)	[7440-47-3]	2,5	%		
	Molybdenum (Mo)	[7439-98-7]	0,22	%		
	Lead (Pb)	[7439-92-1]	0,007	%		
	Sulfur (S)	[7704-34-9]	(<0,01)	%		
	Silicon (Si)	[7440-21-3]	~2	%		
	Tin (Sn)	[7440-31-5]	0,22	%		
	Titanium (Ti)	[7440-32-6]	0,33	%		