

Iron base, Cast Iron, RM, Block 40 mm x 37 mm x 10 mm

Art. ID BAS-LARM 3
Unit block
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,042	%		
	Arsenic (As)	[7440-38-2]	0,092	%		
	Boron (B)	[7440-42-8]	0,003	%		
	Bismuth (Bi)	[7440-69-9]	0,022	%		
	Carbon (C)	[7440-44-0]	~3	%		
	Copper (Cu)	[7440-50-8]	1,2	%		
	Nickel (Ni)	[7440-02-0]	1,8	%		
	Sulfur (S)	[7704-34-9]	(<0,01)	%		
	Silicon (Si)	[7440-21-3]	~2	%		
	Vanadium (V)	[7440-62-2]	0,55	%		