

Austenitic (Ni-RESIST), Iron base, RM, Block 40 mm x 37 mm x 10 mm

Art. ID BAS-NIRM 3
Unit block
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Carbon (C)	[7440-44-0]	2,51	%		
	Cerium (Ce)	[7440-45-1]	0,007	%		
	Chromium (Cr)	[7440-47-3]	2,43	%		
	Copper (Cu)	[7440-50-8]	1	%		
	Manganese (Mn)	[7439-96-5]	0,51	%		
	Niobium (Nb)	[7440-03-1]	0,09	%		
	Nickel (Ni)	[7440-02-0]	17,8	%		
	Phosphorus (P)	[7723-14-0]	0,208	%		
	Sulfur (S)	[7704-34-9]	0,096	%		
	Silicon (Si)	[7440-21-3]	2,21	%		