

**Out of Stock - Item is not available at this time - Density Standard for Medical Computed Tomography**

Art. ID                      NIST-2088  
Unit                         five foam blocks  
Deliverydetails            No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM®) is intended primarily for use in calibrating the density scale of medical computed tomography (CT) machines in the range of the lung. A unit of NIST-2088 consists of five foam blocks held in place by plastic shrink wrap.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Please ask for the current certificate						