

**Out of Stock - Item is not available at this time - Rockwell C Hardness, Low Range**

Art. ID	NIST-2810
Unit	each
Deliverydetails	No Dangerous Good /not restricted

**Description**

This Standard Reference Material (SRM®) is a transfer standard intended for use in the calibration and verification of the performance of Rockwell hardness equipment using the Rockwell C Hardness Scale (HRC). A unit of NIST-2810 is a steel test block, nominally 64 mm in diameter and 15 mm thick, having a polished test surface with a micro-engraved circle 52 mm in diameter. Each SRM unit bears a unique serial number on the edge of the block.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
Please ask for the current certificate						