

Rubidium (Rb) Standard Solution

Art. ID NIST-3145a
Unit 5 x 10 mL
Deliverydetails No Dangerous Good /not restricted

Description

This Standard Reference Material (SRM®) is intended for use as a primary calibration standard for the quantitative determination of rubidium. A unit of NIST-3145a consists of five 10 mL sealed borosilicate glass ampoules of an acidified aqueous solution prepared gravimetrically to contain a known mass fraction of rubidium. The solution contains nitric acid at a volume fraction of approximately 1 %. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Rubidium (Rb)	[7440-17-7]	10,04 ± 0,06	mg/g		