

**Out of Stock - Item is not available at this time - Tungsten (W) Standard Solution**

Art. ID	NIST-3163
Unit	50 mL
Deliverydetails	ADR Limited Quantity (LQ) / AIR Dangerous Good (DG) UN1760 II E2

**Description**

This Standard Reference Material (SRM®) is intended for use as a primary calibration standard for the quantitative determination of tungsten. A unit of NIST-3163 consists of 50 mL of a single element solution in a high-density polyethylene bottle sealed in an aluminized bag. The solution was prepared gravimetrically to contain a known mass fraction of tungsten. The solution contains nitric acid at a concentration of approximately 1.1 mol/L and hydrofluoric acid at a concentration of approximately 1.1 mol/L. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Tungsten (W)	[7440-33-7]	10 ± 0,026	mg/g		