

Creatinine in Frozen Human Urine

Art. ID NIST-3667
Unit 10 mL
Deliverydetails Dry ice shipment /not restricted

Description

This Standard Reference Material (SRM®) is intended primarily for use in evaluating the accuracy of procedures for the determination of creatinine in human urine. It is also intended for use in validating working or secondary reference materials. NIST-3667 was prepared from normal human urine collected from male and female donors, and the creatinine concentration has not been modified. A unit of NIST-3667 consists of one bottle of 10 mL frozen human urine. /// Sample value(s) - please ask for current certificate.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Creatinine	[60-27-5]	613 ± 13	µg/g		
	Creatinine	[60-27-5]	61.8 ± 1.3	mg/dL		