

**Certified Reference Material - Complex Nutrients (Hard Water), effluent (synthetic waste water matrix)**

Art. ID                      ERA-1241  
Unit                            each  
Deliverydetails            ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 III E1

Description

European Union Water Framework Directive Standard. One 15 mL screw cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total Kjeldahl Nitrogen (TKN)		0,5 - 5	mg/L		
	Total nitrogen		1,0 - 20,0	mg/L		
	Total phosphorus		0,5 - 5	mg/L		