

Certified Reference Material - Demand, effluent (synthetic waste water matrix)

Art. ID	ERA-1242
Unit	each
Deliverydetails	ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1789 III E1

Description

European Union Water Framework Directive Standard. One 15 mL screw cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Biochemical Oxygen Demand (BOD)		5,0 - 75,0	mg/L		
	Chemical Oxygen Demand (COD)		10 - 150	mg/L		
	Dissolved Organic Carbon (DOC)		2,0 - 40,0	mg/L		
	Total Organic Carbon (TOC)		2,0 - 40,0	mg/L		