

Certified Reference Material - Solids, effluent (synthetic waste water matrix)

Art. ID ERA-1243
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

European Union Water Framework Directive Standard. One 24 mL screw cap vial yielding 1 liter of solution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total Dissolved Solids at 180 °C (TDS)		10 - 250	mg/L		
	Total suspended solids (TSS)		5,0 - 50	mg/L		