

Certified Reference Material - Solids, effluent (synthetic waste water matrix)

Art. ID ERA-1243
Unit each
Deliverydetails No Dangerous Good /not restricted

Description

European Union Water Framework Directive Standard. One 24 mL screw cap vial yielding 1 liter of solution.

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|------------|---------------------|------|--------|--------|
| | Total Dissolved Solids at 180 °C (TDS) | | 10 - 250 | mg/L | | |
| | Total suspended solids (TSS) | | 5,0 - 50 | mg/L | | |