

Certified Reference Material - Mercury, clean water (synthetic drinking water matrix)

Art. ID ERA-1341
Unit 1 x 15 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN3264 III E1

Description

European Union Water Framework Directive Standard. One 15 mL screw cap vial yields up to 1 liter after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Total mercury		0,1 to 1,2	µg/L		