

**Certified Reference Material - Free Cyanide, clean water (synthetic drinking water matrix)**

Art. ID ERA-1345  
Unit 15 mL vial  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1824 III E1

**Description**

European Union Water Framework Directive Standard. One 15 ml screw cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Free cyanide		5 - 100	µg/L		
	Total cyanide		5 - 100	µg/L		