

**Certified Reference Material - Simple Nutrients (Soft Water), clean water (synthetic drinking water matrix)**

Art. ID                      ERA-1349  
Unit                         each  
Deliverydetails            No Dangerous Good /not restricted

Description

European Union Water Framework Directive Standard. One 15 mL screw cap vial yields up to 2 liters after dilution.

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Ammonium (NH <sub>4</sub> <sup>+</sup> )	[14798-03-9]	0,1 - 1	mg/L		
	Nitrate (NO <sub>3</sub> <sup>-</sup> )		3,0 -60,0	mg/L		
	Nitrite (NO <sub>2</sub> <sup>-</sup> )		0,1 - 1	mg/L		
	Soluble reactive phosphorus (P)		0,5 - 5	mg/L		
	Total oxidised nitrogen (TON)		3,0 -60,0	mg/L		