

**Cast Iron, CRM, High Phosphorus, Block 48 mm x 42 mm x 12 mm**

Art. ID                      BAS-SCRM-661/4  
Unit                         block  
Deliverydetails            No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
|                  | Carbon (C)        | [7440-44-0] | 2,56                | %    |        |        |
|                  | Chromium (Cr)     | [7440-47-3] | ~1                  | %    |        |        |
|                  | Manganese (Mn)    | [7439-96-5] | 0,3                 | %    |        |        |
|                  | Phosphorus (P)    | [7723-14-0] | 0,84                | %    |        |        |
|                  | Sulfur (S)        | [7704-34-9] | 0,068               | %    |        |        |
|                  | Silicon (Si)      | [7440-21-3] | 2,96                | %    |        |        |