

Stainless Steel , Wrought, Certified Reference Material (CRM), disc 38 mm x 15 mm

Art. ID MBH-13X32100A
Unit disc
Deliverydetails No Dangerous Good /not restricted

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|-------------------|-------------|---------------------|------|--------|--------|
| | Aluminium (Al) | [7429-90-5] | 0,0247 | % | | |
| | Boron (B) | [7440-42-8] | 0,0025 | % | | |
| | Carbon (C) | [7440-44-0] | 0,0463 | % | | |
| | Cobalt (Co) | [7440-48-4] | 0,105 | % | | |
| | Chromium (Cr) | [7440-47-3] | 17,39 | % | | |
| | Copper (Cu) | [7440-50-8] | 0,415 | % | | |
| | Manganese (Mn) | [7439-96-5] | 1,52 | % | | |
| | Molybdenum (Mo) | [7439-98-7] | 0,282 | % | | |
| | Nitrogen (N) | [7727-37-9] | 0,0115 | % | | |
| | Niobium (Nb) | [7440-03-1] | 0,0191 | % | | |
| | Nickel (Ni) | [7440-02-0] | 9,32 | % | | |
| | Phosphorus (P) | [7723-14-0] | 0,0298 | % | | |
| | Sulfur (S) | [7704-34-9] | 0,0011 | % | | |
| | Silicon (Si) | [7440-21-3] | 0,498 | % | | |
| | Tin (Sn) | [7440-31-5] | 0,0115 | % | | |
| | Titanium (Ti) | [7440-32-6] | 0,376 | % | | |
| | Vanadium (V) | [7440-62-2] | 0,106 | % | | |
| | Tungsten (W) | [7440-33-7] | 0,021 | % | | |