

**Phenanthrene D10 10 µg/mL in Acetone (Dr. Ehrenstorfer)**

Art. ID DRE-L20920100AC  
Unit 10 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1090 II E2

| Text/Information | Analyte/Parameter | CAS number  | Concentration/Value | Unit  | Method | Source |
|------------------|-------------------|-------------|---------------------|-------|--------|--------|
|                  | Phenanthrene D10  | [1517-22-2] | 10                  | µg/mL |        |        |