

**Quinmerac D4 (methyl(D3)-quinoline-2-D) 100 µg/mL in Acetonitrile**

Art. ID DRE-XA16708100AL  
Unit 1 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Quinmerac D4 (methyl(D3 [n/a] )-quinoline-2-D)		100	µg/mL		