

Internal Standards Mix 33 2000 µg/mL in Toluene

Art. ID DRE-Y08273300TO
Unit 10 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,4-Dichlorobenzene D4	[3855-82-1]	2000	µg/mL		
	Acenaphthene D10	[15067-26-2]	2000	µg/mL		
	Chrysene D12	[1719-03-5]	2000	µg/mL		
	Naphthalene D8	[1146-65-2]	2000	µg/mL		
	Perylene D12	[1520-96-3]	2000	µg/mL		
	Phenanthrene D10	[1517-22-2]	2000	µg/mL		