

PAH-Mix 31 deuterated 1000 µg/mL in Toluene

Art. ID DRE-YA20950031TO
Unit 1 mL
Deliverydetails De Minimis UN1294 II E2

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Acenaphthene D10	[15067-26-2]	1000	µg/mL		
	Chrysene D12	[1719-03-5]	1000	µg/mL		
	Naphthalene D8	[1146-65-2]	1000	µg/mL		
	Perylene D12	[1520-96-3]	1000	µg/mL		
	Phenanthrene D10	[1517-22-2]	1000	µg/mL		