

Silver Steel - Tool Steel Type W-5., Wrought, Certified Reference Material (CRM), disc 40 mm x 15 mm

Art. ID MBH-14X72305A
Unit disc
Deliverydetails No Dangerous Good /not restricted

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Aluminium (Al)	[7429-90-5]	0,0049	%		
	Carbon (C)	[7440-44-0]	1.085	%		
	Chromium (Cr)	[7440-47-3]	0,425	%		
	Copper (Cu)	[7440-50-8]	0,149	%		
	Manganese (Mn)	[7439-96-5]	0,349	%		
	Molybdenum (Mo)	[7439-98-7]	0,0231	%		
	Nitrogen (N)	[7727-37-9]	0,0068	%		
	Nickel (Ni)	[7440-02-0]	0,089	%		
	Phosphorus (P)	[7723-14-0]	0,0128	%		
	Sulfur (S)	[7704-34-9]	0,0028	%		
	Silicon (Si)	[7440-21-3]	0,206	%		
	Tin (Sn)	[7440-31-5]	0,0101	%		
	Vanadium (V)	[7440-62-2]	0,0045	%		