

2-Isopropyl-6-methyl-4-pyrimidinol (methyl-4,5,6-13C4, 99%) 100 µg/mL in Acetonitrile

Art. ID CIL-CLM-4538-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1648 II E2

Description

Chemical purity: 97%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Isopropyl-6-methyl-4-pyrimidinol (methyl-4,5,6-13C4, 99%)	[2483735-79-9]	100	µg/mL		