

2,4,5-Trichlorophenoxyacetic acid (2,4,5-T) (ring-13C6,99%) 100 µg/mL in Methylene chloride

Art. ID CIL-CLM-4551-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1593 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4,5-Trichlorophenoxyacetic acid (2,4,5-T) (ring-13C6,99%)	[1216572-34-7]	100	µg/mL		