

**3-Hydroxyphenanthrene (96% chemical purity) (13C6,98%) 50 µg/mL in Toluene**

Art. ID CIL-CLM-4859-T-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Description

Chemical purity: 96%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	3-hydroxyphenanthrene 9 6% chemical purity(13C6 ,98%)	[1531566-87-6]	50	µg/mL		