

1,2-Dichloroethylene (13C1,99%) (cis/trans mix) (stabilised with hydroquinone) 100 µg/mL in Methanol

Art. ID CIL-CLM-6145-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	1,2-Dichloroethylene (1 3C1,99%) (cis/trans mix)	[2483736-22-5]	100	µg/mL		