

**Cypermethrin (mix of stereoisomers) (phenoxy-13C6,99%) 100 µg/mL in Nonane**

Art. ID CIL-CLM-7292-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cypermethrin (mix of stereoisomers) (phenoxy-13C6,99%)	[52315-07-8]/unlabeled	100	µg/mL		