

trans-Permethrin (phenoxy-13C6,99%) 50 µg/mL in Nonane

Art. ID CIL-CLM-7323-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	trans-Permethrin (phenoxy-13C6,99%)	[51877-74-8]/unlabeled	50	µg/mL		