

2-Hydroxyphenanthrene (13C6,99%) 50 µg/mL in Toluene

Art. ID CIL-CLM-8463-T-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2-Hydroxyphenanthrene (13C6,99%)	[605-55-0]/unlabeled	50	µg/mL		