

Dechlorane Plus® anti-Isomer (13C10,99%) 100 µg/mL in Nonane

Art. ID CIL-CLM-8588-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1920 III E1

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dechlorane Plus anti-Isomer (13C10,99%)		100	µg/mL		