

**Dechlorane Plus® anti (bis-cyclopentene-13C10,99%) 100 µg/mL in Toluene**

Art. ID CIL-CLM-8588-T-1.2  
Unit 1,2 mL  
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1294 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Dechlorane Plus anti (b is-cyclopentene-13C10,9 9%)	[n/a]	100	µg/mL		