

**2,4,6-Trinitrotoluene (15N3,98%,13C7,99%) 1000 µg/mL in Benzene wetted with H2O 33% by weight**

Art. ID CIL-CNLM-3643-S  
Unit 1 mL  
Deliverydetails De Minimis UN1114 II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	2,4,6-Trinitrotoluene ( 15N3,98%,13C7,99%)	[202406-62-0]	1000	µg/mL		