

5,5-Diphenylhydantoin (2-13C,99%, 1,3-15N2,98%) 100 µg/mL in Methanol

Art. ID CIL-CNLM-411-1.2

Unit 1,2 mL

Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) UN1230 (6.1) II E2

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	5,5-Diphenylhydantoin (	[78213-26-0]	100	µg/mL		
	2-13C,99%, 1,3-15N2,98%					
	)					