

Cyanuric acid (U-13C3, 99%, U-15N3,98%+) 100 µg/mL in Water

Art. ID CIL-CNLM-4661-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Chemical purity: 90%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Cyanuric acid (U-13C3,99%, U-15N3,98%+)	[n/a]	100	µg/mL		