

Cyanuric acid (U-13C3, 99%, U-15N3,98%+) 100 µg/mL in Water

Art. ID CIL-CNLM-4661-1.2
Unit 1,2 mL
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ) /not restricted

Description

Chemical purity: 90%

| Text/Information | Analyte/Parameter | CAS number | Concentration/Value | Unit | Method | Source |
|------------------|---|------------|---------------------|-------|--------|--------|
| | Cyanuric acid (U-13C3,99%, U-15N3,98%+) | [n/a] | 100 | µg/mL | | |