

Vinyl chloride (D3,98%) (inhibited with Hydroquinone) (require a Breakseal Flask or Cylinder and Valve at an additional charge)

Art. ID CIL-DLM-167-5
Unit 5 g
Deliverydetails ADR Excepted Quantity (EQ) / AIR Excepted Quantity (EQ)

Description

Chemical purity: 98%

Text/Information	Analyte/Parameter	CAS number	Concentration/Value	Unit	Method	Source
	Vinyl chloride (D3,98%)	[6745-35-3]				